

Use your phone as a microphone

To use your phone as the microphone with MModal Fluency Direct hand held dication.

Go to the App store on your phone and download Solventum Mobile Micorphone



Right click on the gray bar displyaing your name, then click Microphone Setup.





Click the drop down menu next to Microphone and choose Link a mobile device

Microp	phone Setup	? _ >
	Set up your microphone	
	Microphone Select the microphone you want to use for dictation Record button Choose a microphone button or keyboard key to start and stop recording Calibration Calibrate the microphone to your voice and surroundings	Philips® 3500 Philips® 3500 Link a mobile device RECORP
		Advanced Close

You will be give a 6 digit number to enter in the app on your.

up your r	Link a mobile device	? ×		
ophone the microphone you wa	Please download the Fluency Direct Mobile Microphone ap from Google Play or the App Store. Launch the app and enter this 6-digit temporary code when prompted.	p	3500	
rd button e a microphone button	733301	8	CORD	
ration Ite the microphone to y	Waiting for a device			