

Test Name	GW Code	MT Code	Tube Type	Minimum Volume	Special Handling	Misc Notes	Other Valid Container
ABO/RH	WABORH	ABORH	LAV	2 mL peds tube			
ACETAMINOPHEN	6ACTMN	ACTMNA	PST	1 mL			SST, R, GL
ALBUMIN	6ALB	ALBA	PST	1 mL			SST, R, GL
ALCOHOL	6ALCM	ALCMA	PST	1 mL	Site cleanse with Iodine or soap & water		SST, R, GL
ALKALINE PHOSPHATASE	6ALKP	ALKA	PST	1 mL			SST, R, GL
ALT	6ALT	ALTA	PST	1 mL			SST, R, GL
AMMONIA	6AMON	AMONA	LAV	2 mL peds tube	Collect on Ice STAT 15 min to lab		
AMYLASE	6AML	AMLA	PST	Microtainer			SST, R, GL
AST	6AST	ASTA	PST	1 mL			SST, R, GL
A1C	A1CV	A1CV	LAV	Microtainer			
B12 FOLATE	B12FLC	B12FLC	SST/WRAP	2 mL	Protect form Light		R
BETA HYDROBUTYRATE	6BOHB	BOHBA	PST	1 mL			SST, R, GL
BILI ADULT/BABY PANEL	6BILAP	BILAPA	PST	Microtainer	Protect form Light		SST, R, GL
BILIRUBIN DIRECT	6DBIL	DBILA	PST	Microtainer	Protect form Light		SST, R, GL
BILIRUBIN TOTAL	6TBIL	TBILA	PST	Microtainer	Protect form Light		SST, R, GL
BUN	6BUN	BUNA	PST	Microtainer			SST, R, GL
CALCIUM	6CA	CAA	PST	Microtainer			SST, R, GL
CARBAMAZEPINE	CRBAMC	CRBAMC	PST	1 mL		Also known as Tegretol	SST, R, GL
CARBOXYHEMOGLOBIN	6CO	CO	LH GRN	1 mL			PST, Blood Gas
CEA	CEAC	CEAC	SST	1 mL			R
CORTISOL	CORTLC	CORTLC	PST	1 mL			SST, R, GL
CRP	CRPA	CRPA	PST	1 mL			SST, R, GL
CRP-HIGH SENSITIVITY	HSCRPA	HSCRPA	PST	Microtainer			SST, R, GL
CHEM8	6CHEM8	CHEM8A	PST	Microtainer		Also known as BMP	SST, R, GL
QUICK CHEM 8	6CH81 6CH8IP	CHEM8I	LH GRN	Microtainer			
CBC/AUTO DIFF	6CBC	CBCA	LAV	Microtainer			
CKMB	6TMB	TMBC	PST	1 mL			SST, R, GL
COMP	6COMP	COMP	PST	2 Microtainers		Also known as CMP	SST, R, GL
CREATININE	6CRE	CREA	PST	Microtainer			SST, R, GL
CREATININE X-RAY EVALUATION	6CRIR	CRIR	LH GRN	2 mL	STAT to lab in 10 min		PST, Blood Gas
D-DIMER	6DDIME	DDIME	LT BLUE	FILL TO LINE			
DIGOXIN	6DIG	DIGC	PST	1 mL			SST, R, GL
DILANTIN	6PTNC	PTNCC	PST	1 mL			SST, R, GL
ELECTROLYTES	6LYTES	LYTESA	PST	Microtainer			SST, R, GL
FERRITIN	FERC	FERC	PST	1 mL			SST, R, GL
FSH	FSHC	FSHC	PST	1 mL			SST, R, GL
GLUCOSE	6GLU	GLUA	PST	Microtainer			SST, R, GL
GENTAMICIN	6GENTG 6GENPK 6GENRN	GENTGC GENPKC GENRNC	PST	1 mL			SST, R, GL
H&H	6HGB, 6HCT	HGB, HCT	LAV	Microtainer		Included in CBC	
HCG (QUAL & QUANT)	6BHCGC 6BHCGE 6HCGQ	BHCGEC CHCGEC HCGQC	SST	1 mL			R, PST, GL
HEPATITIS ACUTE PANEL	HPAC	HPAC	SST	2 mL			R
HIV	6HIV12	HIV12C	SST	2 mL			R, PST, GL
IGA	IGAA	IGAA	PST	1 mL			SST, R, GL
IGG	IGGA	IGGA	PST	1 mL			SST, R, GL
IGM	IGMA	IGMA	PST	1 mL			SST, R, GL
STAT IONIZED CALCIUM**	6IC	ICA	LH GRN	4 mL	**STAT** STABLE 2HR ROOMTEMP	Whole Blood	PST, Blood Gas

Test Name	GW Code	MT Code	Tube Type	Minimum Volume	Special Handling	Misc Notes	Other Valid Container
ROUTINE (SEND OUT) IONIZED CALCIUM	ICALS	ICALS	SST	4mL	Unopened Spin within one hour. Refrigerate	Requires own tube	
IRON	FEA	FEA	PST	Microtainer		Fasting Preferred	GL
LACTIC ACID	6LA	LACTA	GRAY	Microtainer	Collect on Ice STAT to lab	Must be venous collection	
LUTENIZING HORMONE	LHC	LHC	PST	1 mL			SST, R, GL
LIPASE	6LIP	LIPA	PST	Microtainer			SST, R, GL
LIVER PANEL	6LIVER	LIVERA	PST	Microtainer		Also known as hepatic panel	SST, R, GL
LIPID	6LIPID	LIPIDA	PST	Microtainer		Fasting	SST, R, GL
LITHIUM	LIA	LIA	SST	1 mL			
MAGNESIUM	6MG	MGA	PST	Microtainer			SST, R, GL
METHOTREXATE	6METX	6METX	PST	1 mL			GL
MONO SPOT	6MONO	MONO	LAV	Microtainer			
PHENOBARBITAL	PHNOC	PHNOC	PST	1 mL			SST, R, GL
PFA	6PFA	PFA	LT BLUE	FILL TO LINE (3)	No spin	3 tubes needed	
PHOSPHORUS	6PHOS	PHOA	PST	Microtainer			SST, R, GL
PRO-BNP	6PBNP	PBNPC	PST	Microtainer			SST, R, GL
PROCALCITONIN	6PCAL	PCAL	SST	Microtainer			R, PST, GL
POTASSIUM	6KA	KA	PST	Microtainer			SST, R, GL
PROTIME	6PTIMEN	PTIMEN	LT BLUE	FILL TO LINE			
PSA	PSAC	PSAC	SST	1 mL			R, PST, GL
PTT	6PTT	PTT	LT BLUE	FILL TO LINE			
RETICULOCYTES	6RETC	RETC	LAV	Microtainer			
RENAL PANEL	6RENAL	RENALA	PST	Microtainer			SST, R, GL
RHOGAM	WHRIG	RHIG	LAV	4mL			
SALICYLATE	6SAL	SALA	PST	Microtainer			SST, R, GL
SYPHILIS IGG	SYPHB	SYPHB	SST	1 mL			R
TIBC	TBSATA	TBSATA	PST	1mL			GL
T3	TT3C	TT3C	PST	1 mL			SST, R, GL
T4	T4C	T4C	PST	1 mL			SST, R, GL
TESTOSTERONE	TESTOC	TESTOC	PST	1 mL			SST, R, GL
THEOPHYLLINE	THEOC	THEOC	PST	1 mL			SST, R, GL
THYROID PEROXIDASE ANTIBODY	TPOC	TPOC	SST	1 mL			R, PST, GL
TOBRAMYCIN	TOBTGS TOBPKS TOBRNS	TOBTGS TOBPKS TOBRNS	RED	1mL		Performed at St. Vincent	
TOTAL PROTEIN	6TP	TPA	PST	1 mL			SST, R, GL
TROPONIN	6TNT	TNT	PST	Microtainer			GL
TSH	TSH3C	TSH3C	PST	2 Microtainers			SST, R, GL
TYPE AND CROSSMATCH	WXM	XM	LAV	4 mL		Requires blood bank banding	
TYPE AND SCREEN	WTYSCR	TSYCR	LAV	4 mL		Requires blood bank banding	
URIC ACID	6UIC	URICA	PST	1mL			SST, R, GL
VANCOMYCIN	6VANTG, 6VANPK, 6VANRC	VANTGC VANPKC VANRNC	SST	1mL			R
VALPROIC ACID	6VAL	VALC	PST	1mL			SST, R, GL
VITAMIN D	25HVDC	25HVDC	SST	1mL			
WESTERGREEN SED RATE	6WESR	WESR	LAV	1mL		No capillary collection	

Order of Draw:
1. BLOOD CULTURES
2. LT BLUE
3. RED
4. SST
5. PST/LH GRN
6. LAV
7. GRAY